PROVOX Micron HME

Rx
ONLY

UK
CA

42°C
108°F
RT
2°C
36°F


Figure 1


Figure 2


Figure 3


Figure 4

Prescription Information

CAUTION: United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner. The availability of this product without prescription outside the United States may vary from country to country.

Disclaimer

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

Compatible Products; Kompatible Produkte; Compatibele producten; Produits compatibles ; Prodotti compatibili; Productos compatibles; Produtos compatíveis; Kompatibla produkter; Kompatible produkter; Kompatible produkter; Yhteensopivat tuotteet; Samhæfar vörur; Ûhilduvad tooted; Savietojami izstrādājumi; Suderinami gaminiai; Kompatibilni produkty; Kompatibilis termékek; Kompatibilné produkty; Združljivi izdelki; Zgodne produkty; Produse compatibile; Kompatibilni proizvodi; Kompatibilni proizvodi; Συμβατά προϊόντα; Совместимые изделия; Produk-Produk yang Kompatibel; Uyumlu Ürünler; Համատեղելի արտադրանքներ; Produk yang Serasi; 互換性のある製品; 호환 가능한 제품; 相容產品; 相容产品; ผลิตภัณท์ที่เข้ากันได้; Sản phẩm tương thích; מוצרים תואמים; محصولات سازگار; لمنتجات المتوافق

Table 1

Provox StabiliBase
Provox XtraBase
Provox StabiliBase OptiDerm
Provox FlexiDerm
Provox OptiDerm
Provox LaryTube
Provox LaryButton
Provox HME Cassette Adaptor

Contents

EN - ENGLISH	5
DE - DEUTSCH	8
NL - NEDERLANDS	11
FR - FRANÇAIS	15
IT - ITALIANO	18
ES - ESPAÑOL	21
PT - PORTUGUÊS	25
SV - SVENSKA	28
DA - DANSK	31
NO - NORSK	34
FI - SUOMI	37
IS - ÍSLENSKA	40
ET - EESTI	43
LV - LATVIEŠU	46
LT - LIETUVIŲ KALBA	49
CZ - ČESKY	52
HU - MAGYAR	55
SK - SLOVENČINA	58
SL - SLOVENŠČINA	61
PL - POLSKI	64
RO - ROMÂNĂ	67
HR - HRVATSKI	70
SR - SRPSKI	73
EL - ΕΛΛΗΝΙΚΑ	76
TR - TÜRKÇE	80
HY - ՀԱՅԵՐԵՆԻՍԿ	83
RU - РУССКИЙ	86
ID - BAHASA INDONESIA	90
MS - BAHASA MELAYU	93
JA - 日本語	96
KO - 한국어	99
ZHTW - 繁體中文	102
ZHCN - 简体中文	105
TH - ไทย	108
VI - TIẾNG VIỆT	111
HE - תיבוע	116
AR - العربية	119
FA - فارسی	122

1. Intended Use

Provox Micron HME is a heat and moisture exchanger (HME) and air filtration device for patients breathing through a tracheostoma. Provox Micron HME partially restores lost breathing resistance. For patients with a voice prosthesis or surgical speech fistula, it may also facilitate voicing. Provox Micron HME is intended to be used with the attachment devices in the Provox HME System.

2. CONTRAINDICATIONS

The device should only be used in accordance with the Instructions for Use. Patients without the physical, cognitive, or mental ability required to attach, remove, or operate the device themselves, should not use the device independently and should only use it if they are under sufficient supervision of a clinician or a trained caregiver. The device should not be used by patients with a low tidal volume, as the added dead space may cause CO₂ (carbon dioxide) retention. See Section 3.1 Technical Data.

3. Description of the Device

Provox Micron HME is a Heat and Moisture Exchanger combined with an electrostatic filter. The HME is a foam which contains a salt (calcium chloride). The HME retains the heat and moisture of the exhaled air. When inhaling, the retained heat and moisture in the HME is given back to the lungs. HME use may help improve the function of the lungs and reduce problems with e.g., coughing and mucus production. New users may experience slight discomfort in the beginning related to increased breathing resistance. During the first weeks of use, mucus production may seem to increase. This is normal and means that the mucus is getting thinner and easier to cough up. After a few weeks of HME use, this should stabilize and coughing and mucus production usually decreases.

The Provox Micron HME lid can be pressed down to occlude the stoma in order to speak with a voice prosthesis. When the pressure is released, the lid automatically comes up and the airway passage opens.

Provox Micron HME helps to filter inhaled and exhaled air through consistent normal use. Thereby, small airborne particles, e.g., bacteria, viruses, dust, and pollen are restricted from passing through the device into the lungs upon inspiration (see Technical Data below), while also protecting others nearby during expiration. **Note:** Provox Micron is not intended to be used as a Personal Protective Equipment during work that requires breathing protection.

Information

The products are intended to be stored in the box until removed for use, as the box protects the products and the label on the box contains important information like batch code, use-by date and a unique barcode.

3.1 Technical Data

Pressure drop after 1 h at: (According to ISO 9360)	30 l/min 60 l/min 90 l/min	0.7 hPa 1.9 hPa 3.5 hPa
Pressure drop after 24 h at: (According to ISO 9360)	30 l/min 60 l/min 90 l/min	0.7 hPa 1.8 hPa 3.5 hPa
Moisture output: (According to ISO 9360)	22 mg/L	
Moisture loss at VT=1000 ml: (According to ISO 9360)	22 mg/L	
Bacterial Filtration Efficiency (BFE): (According to ASTM F2101)	≥99%	
Virus Filtration Efficiency (VFE): (According to MIL-M-36954C and ASTM F2101)	≥99%	
HME internal volume (Dead space) (According to ISO 9360)	14.9 ml	

4. WARNING

- Do not exert unintentional physical pressure on the device. Unintentional or accidental closing or compression of the device may obstruct breathing. By putting pressure on the device lid, the airway passage is closed. Inform patient, caregivers, and others of this feature to ensure that they understand the closing feature. Closing the airway to allow voicing is a well-known feature for laryngectomized patient with a voice prosthesis but may be unknown for patients without a voice prosthesis.

5. PRECAUTIONS

- The device provides good protection through consistent normal use as long as there is no air leakage. However, since there are other pathways for e.g., viruses and bacteria to enter the human body, total protection can never be guaranteed.
- The device shall not be used beyond 24 hours, is not reusable, and shall not be rinsed in water or in any other solutions. This will substantially reduce the function of the device and can increase the risk for infections due to bacterial colonization.
- Always test the function of the device prior to use. The lid should immediately return to its open position after releasing it.
- The device should not be used if it impedes directly on the neck, unless a protective layer such as hydrocolloid adhesive or dressing is present, as this may cause skin damage.
- Do not use non-compatible devices or non-original devices as it may cause personal injury or damage to the device. Do not dismantle or alter the device.

- Do not administer medicated nebulizer treatment over the device since the medication can be deposited in the device.
- Do not use humidifiers or heated humidified oxygen via a mask over the tracheostoma while using the device. The device will become too wet. If oxygen therapy is required, use only nonheated humidified oxygen.

6. Instructions for Use

To Attach

Hold the device and press lightly until the device connects to the attachment device. See Fig. 1-2.

To Speak

Press the lid to allow speech. See Fig. 3.

To Remove

Hold the device and pull sideways until the device disconnects. See Fig. 4.

7. Device Lifetime and Disposal

Always follow medical practice and national requirements regarding biohazards when disposing of a used medical device.

- HMEs are single use devices that are designed to provide pulmonary rehabilitation both day and night and are to be replaced in line with daily activities.
- The number and type of HMEs needed during a 24-hour period will vary based on the situations encountered and the frequency of involuntary coughing.
- Typically, one HME is used during the night and about two HMEs during the day. The HME may need to be replaced more often if clogged due to mucus secretions or other contamination. HMEs are not reusable and cannot be rinsed in water or any other solution - this washes out the special salt that is essential for HME effectiveness.

8. Compatible Products

Provox Micron HME is compatible with Provox attachments and accessories. See Table 1.

9. Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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UK Responsible Person
Atos Medical UK Ltd Tottle Road
Cartwright House
Nottingham
Nottinghamshire NG2 1RT
England United Kingdom

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THIS PAGE IS DELIBERATELY LEFT EMPTY

THIS PAGE IS DELIBERATELY LEFT EMPTY

THIS PAGE IS DELIBERATELY LEFT EMPTY

THIS PAGE IS DELIBERATELY LEFT EMPTY

10029, 03
2025-10-30

Atos

atosmedical.com

Atos Medical AB
Kraftgatan 8, SE-242 35 Hörby, Sweden
Tel: +46 (0)415 198 00 • info@atosmedical.com

Atos and the Atos logo are trademarks of Coloplast A/S.
© 2025 Coloplast A/S. All rights reserved.