

# FREEVENT®

XtraCare™  
XtraCare™ *Mini*


Freeevent XtraCare


Freeevent XtraCare Mini


**Figure 1**


**Figure 2**

**Figure 3**

**Figure 4**

## **Prescription information**

**CAUTION:** United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner. The availability of this product without prescription outside the United States may vary from country to country.

## **Disclaimer**

Atos Medical offers no warranty – neither expressed nor implied – to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

## **Patents and Trademarks**

Freevent® is a registered trademark owned by Atos Medical AB, Sweden. XtraCare™ is a trademark of Atos Medical AB. For information about protective rights (e.g. patents), please refer to our web page [www.atosmedical.com/patents](http://www.atosmedical.com/patents).

# Contents

EN - ENGLISH .....	5
DE - DEUTSCH .....	8
NL - NEDERLANDS .....	11
FR - FRANÇAIS .....	14
IT - ITALIANO .....	17
ES - ESPAÑOL .....	20
PT - PORTUGUÊS .....	23
SV - SVENSKA .....	26
DA - DANSK .....	29
NO - NORSK .....	32
FI - SUOMI .....	35
IS - ÍSLENSKA .....	38
ET - EESTI .....	41
LT - LIETUVIŲ KALBA .....	44
CS - ČESKY .....	47
HU - MAGYAR .....	50
SK - SLOVENČINA .....	53
SL - SLOVENŠČINA .....	56
PL - POLSKI .....	59
RO - ROMÂNĂ .....	62
HR - HRVATSKI .....	65
SR - SRPSKI .....	68
EL - ΕΛΛΗΝΙΚΑ .....	71
TR - TÜRKÇE .....	74
HY - ՀԱՅԵՐԵՆ	77
RU - РУССКИЙ .....	80
MS - BAHASA MELAYU .....	83
JA - 日本語 .....	86
KO - 한국어 .....	89
ZHTW - 繁體中文 .....	92
HE - עברית .....	95
AR - العربية .....	98

## 1. Descriptive information

### 1.1 Intended use

Freeevent® XtraCare™ and Freeevent® XtraCare™ Mini are single use Heat and Moisture Exchangers with electrostatic filter (HMEF) that condition and filter inhaled air in patients spontaneously breathing through a tracheostoma.

### 1.2 CONTRAINDICATIONS

Freeevent XtraCare and Freeevent XtraCare Mini are contraindicated for patients who:

- are under any form of mechanical ventilation.
- are unable to handle or remove the device themselves when needed, and who are not under constant supervision of a clinician or a trained caregiver.
- cannot tolerate the added dead space.

Freeevent XtraCare and Freeevent XtraCare Mini are contraindicated if single lumen tracheal tube cannot be reinserted by patient or caregiver after accidentally removed.

### 1.3 Description of the device

The device is available in two versions:

- Freeevent XtraCare.
- Freeevent XtraCare Mini

Both devices are included when referenced as Freeevent XtraCare below.

Freeevent XtraCare is a Heat and Moisture Exchanger combined with an electrostatic filter (HMEF). The HME is impregnated with a hygroscopic salt and conditions the inhaled air. The electrostatic filter reduces the inhalation of particles such as viruses, bacteria, pollen and other particulate matter through the tracheostoma. The electrostatic filter provides effective filtration of small airborne particles, e.g. bacteria, viruses, dust and pollen, irrespective of the direction of airflow through the device. NOTE: Since there are other pathways for pathogens to enter the human body, Freeevent XtraCare can never guarantee total protection. Freeevent XtraCare has a 15 mm ISO connector for connection to a tracheostomy tube. For supplemental oxygen the compatible O<sub>2</sub> adaptor can be connected to Freeevent XtraCare. See section 3.

## 1.4 Technical data

	Freeevent XtraCare	Freeevent XtraCare Mini
Dead space	7.4 ml	6.1 ml
Tidal volume range	>50 ml	30-250 ml
Pressure Drop at*		
- 10 l/min	-	45 Pa
- 30 l/min	70 Pa	150 Pa
- 60 l/min	150 Pa	-
Moisture loss* at		
- VT= 250 ml	-	13.5 mg/l
- VT= 500 ml	18 mg/l	-
- VT= 1000 ml	19.5 mg/l	21.5 mg/l
Filtration efficiency**	At 30 l/min	At 15 l/min
- Virus (VFE)	>99%	>98%
- Bacteria (BFE)	>99%	>99%

\* According to ISO 9360

\*\* Test method adapted from ASTM F2101

## 1.5 WARNINGS

**DO NOT** reuse Freeevent XtraCare on another patient. The device and its accessories are intended for single use only. Reuse in another patient may cause cross contamination.

## 1.6 PRECAUTIONS

- Do not use a damaged or disassembled Freeevent XtraCare as this may cause aspiration of small parts.
- Do not use a soiled or contaminated device as this may cause an infection.
- Do not wash and re-use the HME. Washing the HME impairs the filtering- and HME functions.
- Do not use the device beyond 24 hours after initial use. This can increase the risk for infection due to pathogens.
- Do not administer medicated nebulizer treatment over the device since the medication can become deposited in the device.
- Do not use humidifiers or heated humidified oxygen over the device since the HME will become too wet and may result in an increased work of breathing.

## **2. Instructions for use**

The device is for single use and must be replaced at least every 24 hours, or more often as needed (e.g. if clogged by mucus or wet).

To use Freevent XtraCare, gently attach the device to the tube using a twisting motion (Fig. 1). To remove the device from the tracheostomy tube, hold the tracheostomy tube in place and gently pull the device using a twisting motion.

### **2.1 Disposal**

Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

## **3. Compatible products**

### **Freevent O<sub>2</sub> Adaptor / Freevent O<sub>2</sub> Adaptor Mini:**

Freevent O<sub>2</sub> Adaptor / Freevent O<sub>2</sub> Adaptor Mini are an accessory for patients requiring supplemental oxygen. The adaptor connects oxygen tubing of 1/8 inch (3.2 mm) diameter to the oxygen port.

Attach and detach the O<sub>2</sub> Adaptor to the device as shown in Fig. 2 and Fig. 4. Attach the oxygen tubing to the adaptor as shown in Fig. 3. Make sure the oxygen tube goes beyond the conical edge of the connector (Fig. 3). Flow rates of up to 15 l/min can be used. However, the HME effect will be reduced as the oxygen flow increases. At a typical flow of 3-4 l/min, the HME performance is minimally affected but at 15 l/min HME cannot completely compensate for the extra drying caused by supplementary oxygen.

### **CAUTION:**

- The adaptor is for single use and should be replaced at least every 24 hours, or when it becomes dirty, contaminated or shows any signs of damage.
- Freevent O<sub>2</sub> Adaptor Mini is small and may present a hazard for small children if ingested.

## **4. Reporting**

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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